

CERTIFICATE OF COMPLIANCE

LOT #: 639903

KEY SCIENTIFIC PRODUCTS
1113 E REYNOLDS STREET
STAMFORD, TEXAS 79553
1-800-THE-1-KEY

THIS LOT HAS BEEN TESTED ACCORDING TO THE CRITERIA LISTED
BELOW AND MEETS THE PERFORMANCE STANDARDS AS ESTABLISHED
BY KEY SCIENTIFIC PRODUCTS TO MEET NCCLS GUIDELINES.

PRODUCT: HIPPURATE DISCS 50/BTL
CATALOG #: K710310

LOT #: 639903

PACKAGED DATE = 2024/02/29 EXP DATE = 2026/02/28

CONTROL LOT# 580800

TESTED BY: BLC

	POSITIVE CURRENT	NEGATIVE CURRENT	POSITIVE CONTROL	NEGATIVE CONTROL
12386 19615	+	-	+	-
GROUP B STREP STREPTOCOCCUS PYOGEN				
REAGENT # 614100				
DISK COLOR WHITE INCUBATE FOR 2 HOURS				
POSITIVE: PURPLE COLOR AFTER ADDITION OF DEVELOPER.				
NEGATIVE: NO COLOR CHANGE AFTER ADDITION OF DEVELOPER.				
CHECK HIPPURATE VOL /DISK				
WEIGH 10 BLANK DISKS OF SAME SIZE				
WEIGH 10 FINISHED HIPPURATE DISKS				
DIFFERENCE SHD BE 1.8 TO 2.2 MG PER DISK				