

C E R T I F I C A T E O F C O M P L I A N C E

LOT #: 639905

KEY SCIENTIFIC PRODUCTS
1113 E REYNOLDS STREET
STAMFORD, TEXAS 79553
1-800-THE-1-KEY

THIS LOT HAS BEEN TESTED ACCORDING TO THE CRITERIA LISTED
BELOW AND MEETS THE PERFORMANCE STANDARDS AS ESTABLISHED
BY KEY SCIENTIFIC PRODUCTS TO MEET NCCLS GUIDELINES.

PRODUCT: HIPPURATE DISCS 50/BTL
CATALOG #: K710310

LOT #: 639905

PACKAGED DATE = 2024/08/29 EXP DATE = 2026/08/29

CONTROL LOT# 580800

TESTED BY: BLC

| | | POSITIVE CURRENT | NEGATIVE CURRENT | POSITIVE CONTROL | NEGATIVE CONTROL |
|-------|----------------------|---------------------|---------------------|---------------------|---------------------|
| 12386 | GROUP B STREP | + | | + | |
| 19615 | STREPTOCOCCUS PYOGEN | | - | | - |

REAGENT # 614100

DISK COLOR WHITE INCUBATE FOR 2 HOURS

POSITIVE: PURPLE COLOR AFTER ADDITION OF DEVELOPER.

NEGATIVE: NO COLOR CHANGE AFTER ADDITION OF DEVELOPER.

CHECK HIPPURATE VOL /DISK

WEIGH 10 BLANK DISKS OF SAME SIZE

WEIGH 10 FINISHED HIPPURATE DISKS

DIFFERENCE SHD BE 1.8 TO 2.2 MG PER DISK