

C E R T I F I C A T E   O F   C O M P L I A N C E

LOT #: 580821

KEY SCIENTIFIC PRODUCTS  
1113 E REYNOLDS STREET  
STAMFORD, TEXAS 79553  
1-800-THE-1-KEY

THIS LOT HAS BEEN TESTED ACCORDING TO THE CRITERIA LISTED  
BELOW AND MEETS THE PERFORMANCE STANDARDS AS ESTABLISHED  
BY KEY SCIENTIFIC PRODUCTS TO MEET NCCLS GUIDELINES.

PRODUCT: HIPPURATE DISCS 15/BTL  
CATALOG #: K710315

LOT #: 580821

PACKAGED DATE = 2025/02/20    EXP DATE = 2027/02/20

-----  
CONTROL LOT# 466500

TESTED BY: BLC  
-----

		POSITIVE CURRENT	NEGATIVE CURRENT	POSITIVE CONTROL	NEGATIVE CONTROL
12386	GROUP B STREP	+		+	
19615	STREPTOCOCCUS PYOGEN		-		-

DISK COLOR WHITE                      REAGENT LOT 575100  
POSITIVE: PURPLE COLOR AFTER ADDITION OF DEVELOPER.  
NEGATIVE: NO COLOR CHANGE AFTER ADDITION OF DEVELOPER.  
CHECK HIPPURATE VOL /DISK  
WEIGH 10 BLANK DISKS OF SAME SIZE  
WEIGH 10 FINISHED HIPPURATE DISKS  
DIFFERENCE SHD BE 1.8 TO 2.2 MG PER DISK

RT: 9/2/21 PASSED    3/9/23